

=> file biosis caba caplus embase japio lifesci medline scisearch

=> e bisen prakash/au

```
E1          7      BISEN P S */AU
E2          4      BISEN P S DR/AU
E3          2 --> BISEN PRAKASH/AU
E4         124      BISEN PRAKASH S/AU
E5          1      BISEN PRAKASH S PROF/AU
E6          32      BISEN PRAKASH SINGH/AU
E7          33      BISEN PS/AU
E8           1      BISEN R/AU
E9          24      BISEN R K/AU
E10         12      BISEN R S/AU
E11          1      BISEN RUTH/AU
E12          1      BISEN S/AU
```

=> s e1-e7 and tuberculosis and diagnos?

```
L1          20 ("BISEN P S */AU OR "BISEN P S DR"/AU OR "BISEN PRAKASH"/AU OR
              "BISEN PRAKASH S"/AU OR "BISEN PRAKASH S PROF"/AU OR "BISEN PRAK
              ASH SINGH"/AU OR "BISEN PS"/AU) AND TUBERCULOSIS AND DIAGNOS?
```

=> dup rem l1

PROCESSING COMPLETED FOR L1

```
L2          7 DUP REM L1 (13 DUPLICATES REMOVED)
```

=> s l2 and glycolipid?

```
L3          2 L2 AND GLYCOLIPID?
```

=> d bib ab kwic 1-

YOU HAVE REQUESTED DATA FROM 2 ANSWERS - CONTINUE? Y/(N):y

```
L3  ANSWER 1 OF 2  BIOSIS  COPYRIGHT (c) 2010 The Thomson Corporation  on STN
AN  2005:464785  BIOSIS <<LOGINID::20100824>>
DN  PREV200510248351
TI   ***Glycolipids*** of Mycobacterium ***tuberculosis*** strain H37Rv
    are potential serological markers for ***diagnosis*** of active
    ***tuberculosis*** .
AU  Tiwari, R. P.; Tiwari, Dileep; Garg, Sanjay K.; Chandra, Ramesh;
    ***Bisen, Prakash S.*** [Reprint Author]
CS  Bundelkhand Univ, JC Bose Inst Life Sci, Dept Biotechnol, Jhansi 284218,
    Uttar Pradesh, India
    prakash_bisen@hotmail.com
SO  Clinical and Diagnostic Laboratory Immunology, (MAR 2005) Vol. 12, No. 3,
    pp. 465-473.
    ISSN: 1071-412X.
DT  Article
LA  English
ED  Entered STN: 9 Nov 2005
    Last Updated on STN: 9 Nov 2005
AB  A simple and cost-effective ***diagnostic*** tool (TB Screen Test) for
    the screening of patients with pulmonary and extrapulmonary
    ***tuberculosis*** and for differentiation of those individuals from
    individuals without ***tuberculosis*** , other common infections, and
    healthy controls has been developed. The serological responses of
    purified mycobacterial ***glycolipid*** antigens were examined by a
    liposome agglutination assay. The assay was able to detect very low
    antiglycolipid antibody concentrations in the infected individuals. The
    sera from the ***tuberculosis*** patient group had significantly
```

higher concentrations of antiglycolipid antibody than the sera from uninfected control subjects, with 94% sensitivity and 98.3% specificity.

\*\*\*Glycolipids\*\*\* of Mycobacterium \*\*\*tuberculosis\*\*\* H37Rv antigens were isolated, purified, and characterized. After interchelation with liposome particles, these purified antigens specifically bound to the antiglycolipid antibodies present in the sera of patients with \*\*\*tuberculosis\*\*\*, resulting in the formation of a blue agglutination. This protocol clearly differentiates healthy controls and M. bovis BCG-vaccinated subjects from those with active \*\*\*tuberculosis\*\*\*. The resultant \*\*\*diagnostic\*\*\* tool, the TB Screen Test, is more economical and rapid (4 min) than other currently available products and can be used for the mass screening of a heavily afflicted population.

TI \*\*\*Glycolipids\*\*\* of Mycobacterium \*\*\*tuberculosis\*\*\* strain H37Rv are potential serological markers for \*\*\*diagnosis\*\*\* of active \*\*\*tuberculosis\*\*\*.

AU Tiwari, R. P.; Tiwari, Dileep; Garg, Sanjay K.; Chandra, Ramesh; \*\*\*Bisen, Prakash S.\*\*\* [Reprint Author]

AB A simple and cost-effective \*\*\*diagnostic\*\*\* tool (TB Screen Test) for the screening of patients with pulmonary and extrapulmonary \*\*\*tuberculosis\*\*\* and for differentiation of those individuals from individuals without \*\*\*tuberculosis\*\*\*, other common infections, and healthy controls has been developed. The serological responses of purified mycobacterial \*\*\*glycolipid\*\*\* antigens were examined by a liposome agglutination assay. The assay was able to detect very low antiglycolipid antibody concentrations in the infected individuals. The sera from the \*\*\*tuberculosis\*\*\* patient group had significantly higher concentrations of antiglycolipid antibody than the sera from uninfected control subjects, with 94% sensitivity and 98.3% specificity.

\*\*\*Glycolipids\*\*\* of Mycobacterium \*\*\*tuberculosis\*\*\* H37Rv antigens were isolated, purified, and characterized. After interchelation with liposome particles, these purified antigens specifically bound to the antiglycolipid antibodies present in the sera of patients with \*\*\*tuberculosis\*\*\*, resulting in the formation of a blue agglutination. This protocol clearly differentiates healthy controls and M. bovis BCG-vaccinated subjects from those with active \*\*\*tuberculosis\*\*\*. The resultant \*\*\*diagnostic\*\*\* tool, the TB Screen Test, is more economical and rapid (4 min) than other currently available products and can be. . .

IT Major Concepts  
Infection; Clinical Chemistry (Allied Medical Sciences)

IT Diseases  
\*\*\*tuberculosis\*\*\* : bacterial disease, \*\*\*diagnosis\*\*\*  
\*\*\*Tuberculosis\*\*\* (MeSH)

IT Chemicals & Biochemicals  
\*\*\*glycolipids\*\*\* ; serological markers

IT Methods & Equipment  
TB Screen Test: clinical techniques, \*\*\*diagnostic\*\*\* techniques

ORGN Classifier  
Mycobacteriaceae 08881  
Super Taxa  
Mycobacteria; Actinomycetes and Related Organisms; Eubacteria;  
Bacteria; Microorganisms  
Organism Name  
Mycobacterium bovis (species): pathogen  
Mycobacterium \*\*\*tuberculosis\*\*\* (species): pathogen, strain-H37Rv

Taxa Notes

Bacteria, Eubacteria, Microorganisms

L3 ANSWER 2 OF 2 EMBASE COPYRIGHT (c) 2010 Elsevier B.V. All rights reserved on STN

AN 2007510106 EMBASE <<LOGINID::20100824>>

TI Rapid liposomal agglutination card test for the detection of antigens in patients with active \*\*\*tuberculosis\*\*\* .

AU Tiwari, R.P.

CS Diagnostic Division, Nicholas Piramal India Limited, Pawane, Navi, Mumbai, India.

AU Tiwari, R.P.; Garg, S.K.; \*\*\*Bisen, Prakash S. (correspondence)\*\*\*

CS Institute of Biotechnology and Allied Sciences, Seedling Academy of Design, Technology and Management, Jagatpura, Jaipur, India. psbisen@gmail.com

AU Garg, S.K.

CS Department of Biochemistry, University of Nebraska, Lincoln, NE, United States.

AU Bharmal, R.N.; Kartikeyan, S.

CS Department of Microbiology, Preventive and Social Medicine, Rajiv Gandhi Medical College, Kalwa, Thane, India.

AU \*\*\*Bisen, Prakash S. (correspondence)\*\*\*

CS Bisen Biotech and Biopharma Pvt. Ltd., M-7 Laxmipuram, Transport Nagar, Gwalior 474009, India. psbisen@gmail.com

SO International Journal of Tuberculosis and Lung Disease, (Oct 2007) Vol. 11, No. 10, pp. 1143-1151.

Refs: 30

ISSN: 1027-3719 CODEN: IJTDFO

CY France

DT Journal; Article

FS 015 Chest Diseases, Thoracic Surgery and Tuberculosis

004 Microbiology: Bacteriology, Mycology, Parasitology and Virology

006 Internal Medicine

LA English

SL English; French; Spanish; Castilian

ED Entered STN: 30 Oct 2007

Last Updated on STN: 30 Oct 2007

AB SETTING: A total of 1360 subjects with clinically confirmed pulmonary and extra-pulmonary \*\*\*tuberculosis\*\*\* (TB) and other non-tuberculous conditions. OBJECTIVES: To develop a rapid, sensitive and specific \*\*\*diagnostic\*\*\* test for the detection of the \*\*\*glycolipid\*\*\* antigen of Mycobacterium \*\*\*tuberculosis\*\*\* in a variety of clinical samples. STUDY DESIGN: Affinity-purified rabbit anti- \*\*\*glycolipid\*\*\* antibodies (IgG) were coupled to liposome particles (0.2-0.4 .mu.m) in the presence of 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride and N-hydroxysuccinamide to prepare the working reagent of the TB/M card test. RESULTS: Antibody-conjugated liposomes, when determined with the \*\*\*glycolipid\*\*\* antigens present in the specimens, formed a dark blue agglutination within 4 min. No dumping was observed in samples from normal healthy subjects or patients with other diseases. The test was shown to be effective in detecting \*\*\*glycolipid\*\*\* antigens of M. \*\*\*tuberculosis\*\*\* in clinical samples from patients with active TB

with

as low as 1 ng/ml analytical sensitivity, 97.4% clinical sensitivity and 96.9% specificity. CONCLUSION: The TB/M card test was found to be comparatively economical (4 Indian Rupees or US\$ 0.09/test), rapid (4 min) and seems fairly useful for mass testing of a variety of biological

specimens (cerebrospinal, pleural and synovial fluids, serum, tissue biopsy extract) from patients with tuberculous meningitis, pulmonary TB and other extra-pulmonary TB in endemic countries. .COPYRGT. 2007 The Union.

TI Rapid liposomal agglutination card test for the detection of antigens in patients with active \*\*\*tuberculosis\*\*\* .

AU Tiwari, R.P.; Garg, S.K.; \*\*\*Bisen, Prakash S. (correspondence)\*\*\*

CS Institute of Biotechnology and Allied Sciences, Seedling Academy of Design, Technology and Management, Jagatpura, Jaipur, India

AU \*\*\*Bisen, Prakash S. (correspondence)\*\*\*

CS Bisen Biotech and Biopharma Pvt. Ltd., M-7 Laxmipuram, Transport Nagar, Gwalior 474009, India. psbisen@gmail.com

AB SETTING: A total of 1360 subjects with clinically confirmed pulmonary and extra-pulmonary \*\*\*tuberculosis\*\*\* (TB) and other non-tuberculous conditions. OBJECTIVES: To develop a rapid, sensitive and specific \*\*\*diagnostic\*\*\* test for the detection of the \*\*\*glycolipid\*\*\* antigen of Mycobacterium \*\*\*tuberculosis\*\*\* in a variety of clinical samples. STUDY DESIGN: Affinity-purified rabbit anti- \*\*\*glycolipid\*\*\* antibodies (IgG) were coupled to liposome particles (0.2-0.4 .mu.m) in the presence of 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride and N-hydroxysuccinamide to prepare the working reagent of the TB/M card test. RESULTS: Antibody-conjugated liposomes, when determined with the \*\*\*glycolipid\*\*\* antigens present in the specimens, formed a dark blue agglutination within 4 min. No dumping was observed in samples from normal healthy subjects or patients with other diseases. The test was shown to be effective in detecting \*\*\*glycolipid\*\*\* antigens of M. \*\*\*tuberculosis\*\*\* in clinical samples from patients with active TB

with as low as 1 ng/ml analytical sensitivity, 97.4% clinical sensitivity and.

CT Medical Descriptors:

adolescent

adult

\*agglutination test

\*antigen detection

article

cerebrospinal fluid

controlled study

\*\*\*diagnostic test\*\*\*

\*\*\*extrapulmonary tuberculosis\*\*\*

human

\*\*\*lung tuberculosis\*\*\*

major clinical study

\*\*\*Mycobacterium tuberculosis\*\*\*

pleura fluid

priority journal

school child

sensitivity and specificity

synovial fluid

\*\*\*\*tuberculosis\*\*\*

tuberculous meningitis

1 (3 dimethylaminopropyl) 3 ethylcarbodiimide

amide

antibody conjugate

\*\*\*glycolipid\*\*\*

liposome

n hydroxysuccinamide

tissue extract

=> e tiwari ram/au

E1	1	TIWARI RAKESH VALLABHDAS/AU
E2	28	TIWARI RAKSHA/AU
E3	5 -->	TIWARI RAM/AU
E4	60	TIWARI RAM C/AU
E5	1	TIWARI RAM C DR/AU
E6	15	TIWARI RAM D/AU
E7	11	TIWARI RAM DAS/AU
E8	2	TIWARI RAM KRISHNA/AU
E9	8	TIWARI RAM MOHAN/AU
E10	21	TIWARI RAM P/AU
E11	1	TIWARI RAM P DR/AU
E12	9	TIWARI RAM PARKASH/AU

=> s e1-e12 and tuberculosis and diagnos?

L4 4 ("TIWARI RAKESH VALLABHDAS"/AU OR "TIWARI RAKSHA"/AU OR "TIWARI RAM"/AU OR "TIWARI RAM C"/AU OR "TIWARI RAM C DR"/AU OR "TIWARI RAM D"/AU OR "TIWARI RAM DAS"/AU OR "TIWARI RAM KRISHNA"/AU OR "TIWARI RAM MOHAN"/AU OR "TIWARI RAM P"/AU OR "TIWARI RAM P DR"/AU OR "TIWARI RAM PARKASH"/AU) AND TUBERCULOSIS AND DIAGNOS?

=> dup rem l4

PROCESSING COMPLETED FOR L4

L5 1 DUP REM L4 (3 DUPLICATES REMOVED)

=> d

L5 ANSWER 1 OF 1 BIOSIS COPYRIGHT (c) 2010 The Thomson Corporation on STN  
DUPLICATE 1

AN 2004:37606 BIOSIS <<LOGINID::20100824>>

DN PREV200400038179

TI Analysis of the shotgun expression library of the Mycobacterium  
\*\*\*tuberculosis\*\*\* genome for immunodominant polypeptides: Potential  
use  
in serodiagnosis.

AU Bisen, Prakash S. [Reprint Author]; Garg, Sanjay K.; \*\*\*Tiwari, Ram\*\*\*  
\*\*\* P.\*\*\* ; Tagore, P. Ravindra Nath; Chandra, Ramesh; Karnik, Rucha;  
Thaker,

Nimesh; Desai, Nirav; Ghosh, P. K.; Fraziano, Maurizio; Colizzi, Vittorio  
CS Madhav Institute of Technology and Science, Gwalior, MP, 474 005, India  
prakash\_bisen@hotmail.com

SO Clinical and Diagnostic Laboratory Immunology, (November 2003) Vol. 10,  
No. 6, pp. 1051-1058. print.  
ISSN: 1071-412X (ISSN print).

DT Article

LA English

ED Entered STN: 7 Jan 2004

Last Updated on STN: 7 Jan 2004

=> s l5 and glycolipid?

L6 0 L5 AND GLYCOLIPID?

=> e tiwari ram pramod/au

```

E1          1      TIWARI RAM PRAKAH/AU
E2         20      TIWARI RAM PRAKASH/AU
E3          5 --> TIWARI RAM PRAMOD/AU
E4          1      TIWARI RAMA N/AU
E5          1      TIWARI RAMA NATH/AU
E6          1      TIWARI RAMA SHANKAR/AU
E7          1      TIWARI RAMAPATI/AU
E8          1      TIWARI RAMENDRA K/AU
E9          5      TIWARI RAMESH/AU
E10         5      TIWARI RAMESH CHANDRA/AU
E11         1      TIWARI RAMJI/AU
E12        15      TIWARI RAMMOHAN/AU

```

=> s e3

```
L7          5 "TIWARI RAM PRAMOD"/AU
```

=> dup rem l7

PROCESSING COMPLETED FOR L7

```
L8          1 DUP REM L7 (4 DUPLICATES REMOVED)
```

=> d

```
L8  ANSWER 1 OF 1  BIOSIS  COPYRIGHT (c) 2010 The Thomson Corporation  on STN
    DUPLICATE 1
```

```
AN  2007:385701  BIOSIS <<LOGINID::20100824>>
```

```
DN  PREV200700390349
```

```
TI  Modern approaches to a rapid diagnosis of tuberculosis: Promises and
    challenges ahead.
```

```
AU  ***Tiwari, Ram Pramod*** ; Hattikudur, Narendra S.; Bharmal, Ramesh N.;
    Kartikeyan, S.; Deshmukh, Neeta M.; Bisen, Prakash S. [Reprint Author]
```

```
CS  Seeding Acad Design Technol and Management, Inst Biotechnol and Allied
    Sci, Jaipur 302004, Rajasthan, India
    psbisen@gmail.com
```

```
SO  Tuberculosis (Amsterdam), (MAY 2007) Vol. 87, No. 3, pp. 193-201.
    ISSN: 1472-9792.
```

```
DT  Article
```

```
    General Review; (Literature Review)
```

```
LA  English
```

```
ED  Entered STN: 11 Jul 2007
```

```
    Last Updated on STN: 11 Jul 2007
```

=> s liposom? and tuberculosis and glycolipid and antibod?

```
L9          23 LIPOSOM? AND TUBERCULOSIS AND GLYCOLIPID AND ANTIBOD?
```

=> dup rem l9

PROCESSING COMPLETED FOR L9

```
L10         10 DUP REM L9 (13 DUPLICATES REMOVED)
```

=> s l10 and diagnos?

```
L11         5 L10 AND DIAGNOS?
```

=> d bib ab kwic 1-

YOU HAVE REQUESTED DATA FROM 5 ANSWERS - CONTINUE? Y/(N):y

```
L11  ANSWER 1 OF 5  BIOSIS  COPYRIGHT (c) 2010 The Thomson Corporation  on STN
```

AN 2005:464785 BIOSIS <<LOGINID::20100824>>  
 DN PREV200510248351  
 TI Glycolipids of Mycobacterium \*\*\*tuberculosis\*\*\* strain H37Rv are potential serological markers for \*\*\*diagnosis\*\*\* of active \*\*\*tuberculosis\*\*\* .  
 AU Tiwari, R. P.; Tiwari, Dileep; Garg, Sanjay K.; Chandra, Ramesh; Bisen, Prakash S. [Reprint Author]  
 CS Bundelkhand Univ, JC Bose Inst Life Sci, Dept Biotechnol, Jhansi 284218, Uttar Pradesh, India  
 prakash\_bisen@hotmail.com  
 SO Clinical and Diagnostic Laboratory Immunology, (MAR 2005) Vol. 12, No. 3, pp. 465-473.  
 ISSN: 1071-412X.  
 DT Article  
 LA English  
 ED Entered STN: 9 Nov 2005  
 Last Updated on STN: 9 Nov 2005  
 AB A simple and cost-effective \*\*\*diagnostic\*\*\* tool (TB Screen Test) for the screening of patients with pulmonary and extrapulmonary \*\*\*tuberculosis\*\*\* and for differentiation of those individuals from individuals without \*\*\*tuberculosis\*\*\* , other common infections, and healthy controls has been developed. The serological responses of purified mycobacterial \*\*\*glycolipid\*\*\* antigens were examined by a \*\*\*liposome\*\*\* agglutination assay. The assay was able to detect very low antiglycolipid \*\*\*antibody\*\*\* concentrations in the infected individuals. The sera from the \*\*\*tuberculosis\*\*\* patient group had significantly higher concentrations of antiglycolipid \*\*\*antibody\*\*\* than the sera from uninfected control subjects, with 94% sensitivity and 98.3% specificity. Glycolipids of Mycobacterium \*\*\*tuberculosis\*\*\* H37Rv antigens were isolated, purified, and characterized. After interchelation with \*\*\*liposome\*\*\* particles, these purified antigens specifically bound to the antiglycolipid \*\*\*antibodies\*\*\* present in the sera of patients with \*\*\*tuberculosis\*\*\* , resulting in the formation of a blue agglutination. This protocol clearly differentiates healthy controls and M. bovis BCG-vaccinated subjects from those with active \*\*\*tuberculosis\*\*\* . The resultant \*\*\*diagnostic\*\*\* tool, the TB Screen Test, is more economical and rapid (4 min) than other currently available products and can be used for the mass screening of a heavily afflicted population.  
 TI Glycolipids of Mycobacterium \*\*\*tuberculosis\*\*\* strain H37Rv are potential serological markers for \*\*\*diagnosis\*\*\* of active \*\*\*tuberculosis\*\*\* .  
 AB A simple and cost-effective \*\*\*diagnostic\*\*\* tool (TB Screen Test) for the screening of patients with pulmonary and extrapulmonary \*\*\*tuberculosis\*\*\* and for differentiation of those individuals from individuals without \*\*\*tuberculosis\*\*\* , other common infections, and healthy controls has been developed. The serological responses of purified mycobacterial \*\*\*glycolipid\*\*\* antigens were examined by a \*\*\*liposome\*\*\* agglutination assay. The assay was able to detect very low antiglycolipid \*\*\*antibody\*\*\* concentrations in the infected individuals. The sera from the \*\*\*tuberculosis\*\*\* patient group had significantly higher concentrations of antiglycolipid \*\*\*antibody\*\*\* than the sera from uninfected control subjects, with 94% sensitivity and 98.3% specificity. Glycolipids of Mycobacterium \*\*\*tuberculosis\*\*\* H37Rv antigens were isolated, purified, and characterized. After interchelation with \*\*\*liposome\*\*\* particles, these purified antigens specifically bound to the antiglycolipid \*\*\*antibodies\*\*\* present in

the sera of patients with \*\*\*tuberculosis\*\*\* , resulting in the formation of a blue agglutination. This protocol clearly differentiates healthy controls and M. bovis BCG-vaccinated subjects from those with active \*\*\*tuberculosis\*\*\* . The resultant \*\*\*diagnostic\*\*\* tool, the TB Screen Test, is more economical and rapid (4 min) than other currently available products and can be. . .

IT Major Concepts

Infection; Clinical Chemistry (Allied Medical Sciences)

IT Diseases

\*\*\*tuberculosis\*\*\* : bacterial disease, \*\*\*diagnosis\*\*\*  
 \*\*\*Tuberculosis\*\*\* (MeSH)

IT Chemicals & Biochemicals

glycolipids; serological markers

IT Methods & Equipment

TB Screen Test: clinical techniques, \*\*\*diagnostic\*\*\* techniques

ORGN Classifier

Mycobacteriaceae 08881

Super Taxa

Mycobacteria; Actinomycetes and Related Organisms; Eubacteria;  
 Bacteria; Microorganisms

Organism Name

Mycobacterium bovis (species): pathogen

Mycobacterium \*\*\*tuberculosis\*\*\* (species): pathogen, strain-H37Rv

Taxa Notes

Bacteria, Eubacteria, Microorganisms

L11 ANSWER 2 OF 5 CAPLUS COPYRIGHT 2010 ACS on STN

AN 1997:425363 CAPLUS <<LOGINID::20100824>>

DN 127:32828

OREF 127:6345a,6348a

TI Therapeutic and \*\*\*diagnostic\*\*\* vaccine for the treatment of  
 microbial infections

IN Pascual, David; Bond, Clifford; Burritt, James; Burgess, Don; Glee, Pati;  
 Jutila, John; Jutila, Mark; Bargatze, Robert; Mcfeters, Gordon; Pyle,  
 Barry; Cutler, Jim E.; Han, Yongmoon

PA Research and Development Institute, Inc., USA; Pascual, David; Bond,  
 Clifford; Burritt, James; Burgess, Don; Glee, Pati; Jutila, John; Jutila,  
 Mark; Bargatze, Robert; et al.

SO PCT Int. Appl., 98 pp.

CODEN: PIXXD2

DT Patent

LA English

FAN.CNT 4

	PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
	-----	----	-----	-----	-----
PI	WO 9718790	A2	19970529	WO 1996-US18796	19961121
	WO 9718790	A3	19970731		
	W: AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, IL, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, TJ, TM, TR, TT, UA, UG, US, UZ, VN				
	RW: KE, LS, MW, SD, SZ, UG, AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, ML, MR, NE, SN, TD, TG				
	CA 2238262	A1	19970529	CA 1996-2238262	19961121
	AU 9711226	A	19970611	AU 1997-11226	19961121
	EP 869801	A2	19981014	EP 1996-942049	19961121



EP 869801                      B1            20040121  
R: AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LI, LU, MC, NL,  
PT, SE

JP 2000503630	T	20000328	JP 1997-519932	19961121
AT 258057	T	20040215	AT 1996-942049	19961121
US 20040247611	A1	20041209	US 2004-780650	20040219

PRAI US 1995-7477P            P            19951122  
US 1994-247972                B2            19940523  
US 1995-483558                A2            19950607  
WO 1996-US18796               W            19961121  
US 1998-68935                 B1            19981123

AB Therapeutic peptides, vaccines and \*\*\*diagnostic\*\*\* agents are disclosed for the treatment of pathogenic infections. The agents are capable of binding to mol. address on host cell (e.g. leukocyte, endothelial or epithelial cells, nerve cells), triggering one or more signal transduction pathways and enabling selective pathogen or toxin to traffic through host tissue. The agents are microbial attachment mols. such as adhesive protein, glycoprotein, lectin, carbohydrate, \*\*\*glycolipid\*\*\* .

OSC.G 4            THERE ARE 4 CAPLUS RECORDS THAT CITE THIS RECORD (4 CITINGS)  
RE.CNT 3            THERE ARE 3 CITED REFERENCES AVAILABLE FOR THIS RECORD  
ALL CITATIONS AVAILABLE IN THE RE FORMAT

TI Therapeutic and \*\*\*diagnostic\*\*\* vaccine for the treatment of microbial infections

AB Therapeutic peptides, vaccines and \*\*\*diagnostic\*\*\* agents are disclosed for the treatment of pathogenic infections. The agents are capable of binding to mol. address on host. . . or toxin to traffic through host tissue. The agents are microbial attachment mols. such as adhesive protein, glycoprotein, lectin, carbohydrate, \*\*\*glycolipid\*\*\* .

ST microbial adhesion mol vaccine \*\*\*diagnostic\*\*\* ; monoclonal \*\*\*antibody\*\*\* microbial antigen therapeutic \*\*\*diagnostic\*\*\*

IT Agglutinins and Lectins  
RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(C-type (calcium-dependent type); vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(CD41a; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(CD49; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Gene, animal  
RL: BSU (Biological study, unclassified); BIOL (Biological study)  
(Cdc42, CTP-binding protein; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Selectins  
RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(E-; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Proteins, specific or class

RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(GTP-binding; vaccine comprising microbial adhesion mol. antigen as  
therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Cell adhesion molecules  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(ICAM-1 (intercellular adhesion mol. 1); vaccine comprising microbial  
adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for  
microbial infections)

IT Cell adhesion molecules  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(ICAM-2 (intercellular adhesion mol. 2); vaccine comprising microbial  
adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for  
microbial infections)

IT Cell adhesion molecules  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(ICAM-3 (intercellular adhesion mol. 3); vaccine comprising microbial  
adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for  
microbial infections)

IT Selectins  
RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL  
(Biological study); USES (Uses)  
(L-; vaccine comprising microbial adhesion mol. antigen as therapeutic  
and \*\*\*diagnostic\*\*\* for microbial infections)

IT Cell adhesion molecules  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(MAdCAM-1; vaccine comprising microbial adhesion mol. antigen as  
therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Cell adhesion molecules  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(N-CAM; vaccine comprising microbial adhesion mol. antigen as  
therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Selectins  
RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL  
(Biological study); USES (Uses)  
(P-; vaccine comprising microbial adhesion mol. antigen as therapeutic  
and \*\*\*diagnostic\*\*\* for microbial infections)

IT Cell adhesion molecules  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(PECAM-1; vaccine comprising microbial adhesion mol. antigen as  
therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Toxins  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(Shiga-like toxin, Escherichia coli; vaccine comprising microbial  
adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for  
microbial infections)

IT Antigens  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(Streptococcal SA I/II; vaccine comprising microbial adhesion mol.  
antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial  
infections)

IT Cell adhesion molecules  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(VAM-1; vaccine comprising microbial adhesion mol. antigen as  
therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL  
(Biological study); USES (Uses)

(VLA; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Intestine  
(adhesion of Escherichia coli; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT \*\*\*Diagnosis\*\*\*  
(agents; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Macrophage  
(alveolar; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Lung  
(alveolus, macrophage; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(antigens Mac-1 (macrophage 1); vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Epithelium  
Respiratory tract  
(cells; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Peptides, biological studies  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(conjugates; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Blood vessel  
(endothelium, cells; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Toxins  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(enterotoxins, Escherichia coli; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Clostridium botulinum  
Clostridium tetani  
(exotoxin; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Toxins  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(exotoxins; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT G proteins (guanine nucleotide-binding proteins)  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(gene TC4; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT G proteins (guanine nucleotide-binding proteins)  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(gene rab; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT G proteins (guanine nucleotide-binding proteins)  
RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
(gene rac; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (gp150.95; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Digestive tract  
 Digestive tract  
 (hemorrhage, microorganism causing; vaccine comprising microbial  
 adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for  
 microbial infections)

IT Urogenital tract  
 (infection; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL  
 (Biological study); USES (Uses)  
 (leucam; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Lung  
 (macrophage; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Sialic acids  
 RL: BSU (Biological study, unclassified); BIOL (Biological study)  
 (microbial adhesion mol. contg.; vaccine comprising microbial adhesion  
 mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial  
 infections)

IT Glycopeptides  
 Peptides, biological studies  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (microbial adhesion mol.; vaccine comprising microbial adhesion mol.  
 antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial  
 infections)

IT Agglutinins and Lectins  
 Ligands  
 RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL  
 (Biological study); USES (Uses)  
 (microbial; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT \*\*\*Antibodies\*\*\*  
 RL: BPN (Biosynthetic preparation); THU (Therapeutic use); BIOL  
 (Biological study); PREP (Preparation); USES (Uses)  
 (monoclonal; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Pharynx  
 (nasopharynx, epithelium and endothelium; vaccine comprising microbial  
 adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for  
 microbial infections)

IT Nerve  
 (neuron; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Infection  
 (nosocomial; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Peptides, biological studies  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (oligopeptides; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Urinary tract

(pathogen; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Fungi  
 (phycomycetous; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Paramagnetic materials  
 (superparamagnetic, beads; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Escherichia coli  
 (uropathjogenic; vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Animal cell  
 Animal tissue  
 Aspergillus  
 B cell (lymphocyte)  
 Bacteriophage  
 Blastomyces  
 Bordetella pertussis  
 Brucella  
 Candida  
 Candida albicans  
 Cell adhesion  
 Chlamydia  
 Coccidioides  
 Coliphage M13  
 Corynebacterium diphtheriae  
 Cowpea mosaic virus  
 Cryptococcus (fungus)  
 Cryptosporidium  
 \*\*\*Diagnosis\*\*\*  
 Entamoeba histolytica  
 Enterobacter aerogenes  
 Eukaryote (Eukaryotae)  
 Francisella tularensis  
 Fungi  
 Genetic vectors  
 Giardia lamblia  
 Haemophilus influenzae  
 Hafnia alvei  
 Hantavirus  
 Helicobacter pylori  
 Hepatitis virus  
 Histoplasma  
 Human adenovirus  
 Human coxsackievirus  
 Human herpesvirus  
 Human immunodeficiency virus  
 Human poliovirus  
 Influenza A virus  
 Influenza B virus  
 Influenza C virus  
 Klebsiella pneumoniae  
 Legionella  
 Leishmania  
 Leukocyte  
 \*\*\*Liposomes\*\*\*

Measles virus  
 Microorganism  
 Mumps virus  
 Mycobacterium \*\*\*tuberculosis\*\*\*  
 Mycoplasma pneumoniae  
 Neisseria gonorrhoeae  
 Neisseria meningitidis  
 Organ, animal  
 Parasite  
 Pathogen  
 Pilus  
 Plasmodium berghei  
 Plasmodium falciparum  
 Prokaryote  
 Proteus (bacterium)  
 Pseudomonas  
 Pseudomonas aeruginosa  
 Rhinovirus  
 Rubella virus  
 Salmonella  
 Salmonella typhi  
 Salmonella typhimurium  
 Shigella  
 Signal transduction, biological  
 Staphylococcus  
 Streptococcus  
 T cell (lymphocyte)  
 Treponema pallidum  
 Trichomonas vaginalis  
 Tritrichomonas foetus  
 Trypanosoma  
 Vaccines  
 Vibrio cholerae  
 Yersinia enterocolitica  
 Yersinia pestis  
 Yersinia pseudotuberculosis  
     (vaccine comprising microbial adhesion mol. antigen as therapeutic and  
     \*\*\*diagnostic\*\*\* for microbial infections)  
 IT Cell adhesion molecules  
     RL: ADV (Adverse effect, including toxicity); BSU (Biological study,  
     unclassified); THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
     (vaccine comprising microbial adhesion mol. antigen as therapeutic and  
     \*\*\*diagnostic\*\*\* for microbial infections)  
 IT Antigens  
     Carbohydrates, biological studies  
     Glycolipids  
     Glycoproteins, general, biological studies  
     Integrins  
     Selectins  
     RL: BSU (Biological study, unclassified); THU (Therapeutic use); BIOL  
     (Biological study); USES (Uses)  
     (vaccine comprising microbial adhesion mol. antigen as therapeutic and  
     \*\*\*diagnostic\*\*\* for microbial infections)  
 IT ADP ribosylation factor  
     Adhesins  
     Antiserums  
     Chemokines

Cytokines  
 Glycoconjugates  
 Immunoglobulins  
 LFA-1 (antigen)  
 Ras proteins  
 Rho protein (G protein)  
 Toxins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (vaccine comprising microbial adhesion mol. antigen as therapeutic and  
 \*\*\*diagnostic\*\*\* for microbial infections)

IT Cell adhesion molecules  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (vascular or VCAM; vaccine comprising microbial adhesion mol. antigen  
 as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (.alpha.v; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (.alpha.1.beta.1; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (.alpha.2.beta.1; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (.alpha.3.beta.1; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (.alpha.4.beta.1; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (.alpha.5.beta.1; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT Integrins  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (.alpha.6.beta.1; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT 72146-52-2, Mutan  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (Streptococcal; vaccine comprising microbial adhesion mol. antigen as  
 therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT 9005-32-7, Alginic acid  
 RL: BUU (Biological use, unclassified); BIOL (Biological study); USES  
 (Uses)  
 (gel; vaccine comprising microbial adhesion mol. antigen as therapeutic  
 and \*\*\*diagnostic\*\*\* for microbial infections)

IT 59-23-4, Galactose, biological studies 63-42-3, Lactose 131-48-6,  
 N-Acetylneuraminic acid 1811-31-0, N-Acetylgalactosamine 2438-80-4,  
 Fucose 3416-24-8, Glucosamine 3458-28-4, Mannose 7512-17-6  
 7535-00-4, Galactosamine  
 RL: BSU (Biological study, unclassified); BIOL (Biological study)  
 (microbial adhesion mol. contg.; vaccine comprising microbial adhesion

mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

IT 29350-58-1, PNAd-1  
 RL: THU (Therapeutic use); BIOL (Biological study); USES (Uses)  
 (vaccine comprising microbial adhesion mol. antigen as therapeutic and \*\*\*diagnostic\*\*\* for microbial infections)

L11 ANSWER 3 OF 5 EMBASE COPYRIGHT (c) 2010 Elsevier B.V. All rights reserved on STN

AN 2007510106 EMBASE <<LOGINID::20100824>>

TI Rapid \*\*\*liposomal\*\*\* agglutination card test for the detection of antigens in patients with active \*\*\*tuberculosis\*\*\* .

AU Tiwari, R.P.

CS Diagnostic Division, Nicholas Piramal India Limited, Pawane, Navi, Mumbai, India.

AU Tiwari, R.P.; Garg, S.K.; Bisen, Prakash S. (correspondence)

CS Institute of Biotechnology and Allied Sciences, Seedling Academy of Design, Technology and Management, Jagatpura, Jaipur, India. psbisen@gmail.com

AU Garg, S.K.

CS Department of Biochemistry, University of Nebraska, Lincoln, NE, United States.

AU Bharmal, R.N.; Kartikeyan, S.

CS Department of Microbiology, Preventive and Social Medicine, Rajiv Gandhi Medical College, Kalwa, Thane, India.

AU Bisen, Prakash S. (correspondence)

CS Bisen Biotech and Biopharma Pvt. Ltd., M-7 Laxmipuram, Transport Nagar, Gwalior 474009, India. psbisen@gmail.com

SO International Journal of Tuberculosis and Lung Disease, (Oct 2007) Vol. 11, No. 10, pp. 1143-1151.

Refs: 30  
 ISSN: 1027-3719 CODEN: IJTDFO

CY France

DT Journal; Article

FS 015 Chest Diseases, Thoracic Surgery and Tuberculosis  
 004 Microbiology: Bacteriology, Mycology, Parasitology and Virology  
 006 Internal Medicine

LA English

SL English; French; Spanish; Castilian

ED Entered STN: 30 Oct 2007  
 Last Updated on STN: 30 Oct 2007

AB SETTING: A total of 1360 subjects with clinically confirmed pulmonary and extra-pulmonary \*\*\*tuberculosis\*\*\* (TB) and other non-tuberculous conditions. OBJECTIVES: To develop a rapid, sensitive and specific \*\*\*diagnostic\*\*\* test for the detection of the \*\*\*glycolipid\*\*\* antigen of Mycobacterium \*\*\*tuberculosis\*\*\* in a variety of clinical samples. STUDY DESIGN: Affinity-purified rabbit anti- \*\*\*glycolipid\*\*\* \*\*\*antibodies\*\*\* (IgG) were coupled to \*\*\*liposome\*\*\* particles (0.2-0.4 .mu.m) in the presence of 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride and N-hydroxysuccinamide to prepare the working reagent of the TB/M card test. RESULTS: \*\*\*Antibody\*\*\* -conjugated \*\*\*liposomes\*\*\* , when determined with the \*\*\*glycolipid\*\*\* antigens present in the specimens, formed a dark blue agglutination within 4 min. No dumping was observed in samples from normal healthy subjects or patients with other diseases. The test was shown to be effective in detecting \*\*\*glycolipid\*\*\* antigens of M. \*\*\*tuberculosis\*\*\* in clinical samples from patients with active TB with as low as 1 ng/ml



analytical sensitivity, 97.4% clinical sensitivity and 96.9% specificity.  
CONCLUSION: The TB/M card test was found to be comparatively economical (4 Indian Rupees or US\$ 0.09/test), rapid (4 min) and seems fairly useful for mass testing of a variety of biological specimens (cerebrospinal, pleural and synovial fluids, serum, tissue biopsy extract) from patients with tuberculous meningitis, pulmonary TB and other extra-pulmonary TB in endemic countries. .COPYRGT. 2007 The Union.

TI Rapid \*\*\*liposomal\*\*\* agglutination card test for the detection of antigens in patients with active \*\*\*tuberculosis\*\*\* .

AB SETTING: A total of 1360 subjects with clinically confirmed pulmonary and extra-pulmonary \*\*\*tuberculosis\*\*\* (TB) and other non-tuberculous conditions. OBJECTIVES: To develop a rapid, sensitive and specific \*\*\*diagnostic\*\*\* test for the detection of the \*\*\*glycolipid\*\*\* antigen of Mycobacterium \*\*\*tuberculosis\*\*\* in a variety of clinical samples. STUDY DESIGN: Affinity-purified rabbit anti- \*\*\*glycolipid\*\*\* \*\*\*antibodies\*\*\* (IgG) were coupled to \*\*\*liposome\*\*\* particles (0.2-0.4 .mu.m) in the presence of 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride and N-hydroxysuccinamide to prepare the working reagent of the TB/M card test. RESULTS: \*\*\*Antibody\*\*\* -conjugated \*\*\*liposomes\*\*\* , when determined with the \*\*\*glycolipid\*\*\* antigens present in the specimens, formed a dark blue agglutination within 4 min. No dumping was observed in samples from normal healthy subjects or patients with other diseases. The test was shown to be effective in detecting \*\*\*glycolipid\*\*\* antigens of M. \*\*\*tuberculosis\*\*\* in clinical samples from patients with active TB with as low as 1 ng/ml analytical sensitivity, 97.4% clinical sensitivity and. . .

CT Medical Descriptors:  
adolescent  
adult  
\*agglutination test  
\*antigen detection  
article  
cerebrospinal fluid  
controlled study  
\*\*\*diagnostic test\*\*\*  
\*\*\*extrapulmonary tuberculosis\*\*\*  
human  
\*\*\*lung tuberculosis\*\*\*  
major clinical study  
\*\*\*Mycobacterium tuberculosis\*\*\*  
pleura fluid  
priority journal  
school child  
sensitivity and specificity  
synovial fluid  
\*\*\*\*tuberculosis\*\*\*  
tuberculous meningitis  
1 (3 dimethylaminopropyl) 3 ethylcarbodiimide  
amide  
\*\*\*antibody conjugate\*\*\*  
\*\*\*glycolipid\*\*\*  
\*\*\*liposome\*\*\*  
n hydroxysuccinamide  
tissue extract

AN 0006205491 EMBASE <<LOGINID::20100824>>  
 CP MEDLINE.RTM. is the source for the citation and abstract of this record.  
 TI [Major trends in lipid immunochemistry].  
 Osnovnye napravleniia immunokhimii lipidov..  
 AU Shvets, V.I. (correspondence); Krasnopol'skii, I.M.  
 SO Ukrainskii biokhimicheskii zhurnal, (1984 May-Jun) Vol. 56, No. 3, pp.  
 254-263.  
 ISSN: 0201-8470  
 CY Russian Federation  
 DT Journal; Article  
 FS MEDLINE  
 LA Russian  
 ED Entered STN: Mar 2010  
 Last Updated on STN: Mar 2010  
 AB Data are presented on immunochemical properties of lipids, the most  
 important group of biologically active substances. Problems on antigenic,  
 immunogenic and adjuvant activities of lipids are considered. A possible  
 use of lipid antigens for \*\*\*diagnosis\*\*\* of different infectious  
 diseases is demonstrated and main principles of their construction are  
 suggested. Data are available on immunogenicity of phospho- and  
 \*\*\*glycolipid\*\*\* mixtures as well as on practical application of the  
 obtained \*\*\*antibodies\*\*\*. Guidelines for the use of immunochemical  
 properties of lipids are outlined.  
 AB . . . active substances. Problems on antigenic, immunogenic and  
 adjuvant activities of lipids are considered. A possible use of lipid  
 antigens for \*\*\*diagnosis\*\*\* of different infectious diseases is  
 demonstrated and main principles of their construction are suggested.  
 Data are available on immunogenicity of phospho- and \*\*\*glycolipid\*\*\*  
 mixtures as well as on practical application of the obtained  
 \*\*\*antibodies\*\*\*. Guidelines for the use of immunochemical properties  
 of lipids are outlined.  
 CT Medical Descriptors:  
 animal  
 article  
 brain  
 cattle  
 human  
 immunization  
 immunology  
 \*\*\*lung tuberculosis: DI, diagnosis\*\*\*  
 \*\*\*schistosomiasis: DI, diagnosis\*\*\*  
 serology  
 syphilis serology  
 cardiolipin  
 \*\*\*diagnostic agent\*\*\*  
 \*epitope: AN, drug analysis  
 immunological adjuvant: AD, drug administration  
 \*lipid: AD, drug administration  
 \*\*\*liposome: AD, drug administration\*\*\*  
 L11 ANSWER 5 OF 5 SCISEARCH COPYRIGHT (c) 2010 The Thomson Corporation on  
 STN  
 AN 2008:599597 SCISEARCH <<LOGINID::20100824>>  
 GA The Genuine Article (R) Number: 292ZK  
 TI A novel application of affinity biosensor technology to detect  
 \*\*\*antibodies\*\*\* to mycolic acid in \*\*\*tuberculosis\*\*\* patients  
 AU Verschoor, Jan A. (Reprint)

CS Univ Pretoria, Dept Biochem, ZA-0002 Pretoria, South Africa (Reprint)  
AU Thanyani, Simon T.; Roberts, Vanessa; Siko, D. Gilbert R.; Vrey, Pieter  
CS E-mail: jan.verschoor@up.ac.za  
CYA South Africa  
SO JOURNAL OF IMMUNOLOGICAL METHODS, (20 MAR 2008) Vol. 332, No. 1-2, pp.  
61-72.

ISSN: 0022-1759.

PB ELSEVIER SCIENCE BV, PO BOX 211, 1000 AE AMSTERDAM, NETHERLANDS.

DT Article; Journal

LA English

REC Reference Count: 44

ED Entered STN: 15 May 2008

Last Updated on STN: 15 May 2008

\*ABSTRACT IS AVAILABLE IN THE ALL AND IALL FORMATS\*

AB \*\*\*Tuberculosis\*\*\* has re-emerged as a global health problem due to co-infection with HIV and the emergence of drug-resistant strains of Mycobacterium \*\*\*tuberculosis\*\*\*. HIV co-infection introduced a 30% underestimation in TB \*\*\*diagnosis\*\*\* based on sputum analysis, calling for a reliable and fast serodiagnostic assay to assist in the management of TB in HIV-burdened populations. Serodiagnosis with mycobacterial lipid cell wall antigens gave promising results, in particular with LAM and cord factor. Free mycolic acids have also been considered because they are unique in structure to each species of Mycobacterium and can be economically extracted and purified. In a standard immunoassay such as ELISA, however, an unacceptable number of false positive and false negative test results were obtained. Here we report a much improved biosensor method to detect \*\*\*antibodies\*\*\* to mycolic acids in patient serum as surrogate markers of active \*\*\*tuberculosis\*\*\*. Mycolic acid (MA) \*\*\*liposomes\*\*\* were immobilized on a non-derivatized twin-celled biosensor cuvette and blocked with saponin. A high dilution of serum was used to calibrate the binding signal of the two cells, followed by contact with patient serum at a lesser dilution, but pre-incubated with either antigen-carrying, or empty \*\*\*liposomes\*\*\*. The serum, or the protein A purified IgG thereof,

from

sputum-positive \*\*\*tuberculosis\*\*\* patients could be inhibited from binding to the MA in the biosensor by prior incubation with MA-containing \*\*\*liposomes\*\*\*. The accuracy of the inhibition test was 84% if HIV-positive patients for whom a negative TB sputum analyses could not be relied upon to serve as a reference standard were excluded. If biosensor technology could be made suitable for high throughput screening, then it may provide the solution to the serodiagnosis of \*\*\*tuberculosis\*\*\* against a background of HIV (C) 2007 Elsevier B.V. All rights reserved.

TI A novel application of affinity biosensor technology to detect

\*\*\*antibodies\*\*\* to mycolic acid in \*\*\*tuberculosis\*\*\* patients

AB \*\*\*Tuberculosis\*\*\* has re-emerged as a global health problem due to co-infection with HIV and the emergence of drug-resistant strains of Mycobacterium \*\*\*tuberculosis\*\*\*. HIV co-infection introduced a 30% underestimation in TB \*\*\*diagnosis\*\*\* based on sputum analysis, calling for a reliable and fast serodiagnostic assay to assist in the management of TB in. . . of false positive and false negative test results were obtained. Here we report a much improved biosensor method to detect \*\*\*antibodies\*\*\* to mycolic acids in patient serum as surrogate markers of active \*\*\*tuberculosis\*\*\*. Mycolic acid (MA) \*\*\*liposomes\*\*\* were immobilized on a non-derivatized twin-celled biosensor cuvette and blocked with saponin. A high dilution of serum was used to. . . the two cells, followed by contact with patient serum at a

lesser dilution, but pre-incubated with either antigen-carrying, or empty  
\*\*\*liposomes\*\*\* . The serum, or the protein A purified IgG thereof,  
from sputum-positive \*\*\*tuberculosis\*\*\* patients could be inhibited from  
binding to the MA in the biosensor by prior incubation with MA-containing  
\*\*\*liposomes\*\*\* . The accuracy of the inhibition test was 84% if  
HIV-positive patients for whom a negative TB sputum analyses could not. .  
. biosensor technology could be made suitable for high throughput  
screening, then it may provide the solution to the serodiagnosis of  
\*\*\*tuberculosis\*\*\* against a background of HIV (C) 2007 Elsevier B.V.  
All rights reserved.

ST Author Keywords: \*\*\*antibodies\*\*\* ; mycolic acids; biosensor;  
Mycobacterium \*\*\*tuberculosis\*\*\* ; serodiagnosis

STP KeyWords Plus (R): MYCOBACTERIUM- \*\*\*TUBERCULOSIS\*\*\* ; PULMONARY  
\*\*\*TUBERCULOSIS\*\*\* ; IMMUNE-RESPONSES; PROTEIN-A; BIOMOLECULAR  
INTERACTIONS; SEROLOGICAL \*\*\*DIAGNOSIS\*\*\* ; \*\*\*GLYCOLIPID\*\*\*  
ANTIGEN; OPTICAL BIOSENSOR; RESONANT MIRROR; IGG \*\*\*ANTIBODY\*\*\*